





Technical Data Sheet

2XYT Agar RDM-2XYT-01

Principle

The tryptone serves the source of nitrogen, yeast extract serves the source of nitrogen and vitamins, and sodium chloride balances the osmotic balance. Agar is solidifying agent. The 2xYT agar supports the vigorous growth of $E.\ coli$, necessary for the cultivation of microphages.

Use: Recommended for maintenance and cultivating recombinant strains of *Escherichia coli* and propagation of bacteriophage.

Contents*

Ingredients	Gram/Litre
Tryptone	16.000
Yeast Extract	10.000
Sodium Chloride	5.000
Agar	15.000
pH at 25°C	7.0 ± 0.2

^{*} Formula adjusted for optimum performance and parameters

Directions: Dissolve 46.00 grams in 1000 ml distilled water. Boil to dissolve the medium completely and sterilize by autoclaving at 15 lbs. pressure (121°C) for 15 min, cool it to 42-45 °C and distribute aseptically. Ensure complete solidification and inoculate test sample aseptically.

Specimens types analyzed

Pharmaceutical samples, clinical and non-clinical samples etc.

Precautions to be taken

These microbial media are intended for the in-vitro use only. All the handling, experiments, storage, and discarding should be performed with the help of skilled and knowledgeable technicians and as per the established guidelines. The material should be disposed only after proper sterilization by autoclaving. Please go through the MSDS of the media to avoid any accidents or in emergency.

Performance and Evaluation

The expected performance of the medium is liable to use as per the direction on the label when stored at optimum conditions and within expiry date.

Quality Control

Quanty Control		
Appearance	Beige colored free flowing, homogeneous powder	
Reaction of 4.6% solution	7.00 ±0.2 at 25 °C	
pH	6.80- 7.20	
Gelling strength	Firm comparable with 1.5% agar gel	
Color and clarity of ready medium	Cream to yellow colored opalescent gel	
Growth Promotion properties	Best at ≤ 100 CFU at 32-37 °C for 18-72 h	
Indicative properties	Optimum at ≤ 100 CFU at 32-37 °C for 18-48 h	
Negative control	Performed using sterile distilled water	

Different Microbial Response

Organism	ATCC	Inoculum	Growth	Recovery %	Incubation
Escherichia coli	8739	50-100	Luxurious	70-80 %	33-37 °C, 18-48 h

Storage and Shelf Life

Hygroscopic; keep container tightly closed. Store in cool dry place.

Disposal: To avoid the contamination or propagation of any hazardous microbes the used, unusable or modified preparation of this product must be disposed after autoclaving after completion of task.

Reference

- 1. Atlas, R. M. (2005). Handbook of media for environmental microbiology. CRC press.
- 2. Davis, L.G., M.D. Dibner and J. F. Battey, (1986), *Basic methods in molecular biology*, Elsevier, New York.
- 3. *Difco Manual* (1998). 11th Edition. Difco Laboratories., Division of Becton Dickinson and Company, Sparks, Maryland, USA.
- 4. Rand, M. C., Arnold E. Greenberg, and Michael J. Taras, (1976), *Standard methods for the examination of water and wastewater*. Prepared and published jointly by American Public Health Association, American Water Works Association, and Water Pollution Control Federation.

Disclaimer

The information contained in the technical data sheet is to the best of our knowledge is accurate and true based on the research and development work carried out by <code>ReadyMED</code>[®], Chaitanya Agro Biotech, Malkapur, Maharashtra. The products are neither intended for any therapeutic use for animal or human nor for any other <code>in-vivo</code> applications. The <code>ReadyMED</code>[®] products are only meant to be used for the laboratory, diagnostic, research, or further manufacturing purpose only. These technical outcomes should not be considered as the warranty of any kind expressed or implied, and no liability is accepted for infringement of any patent.

CHAITANYA AGRO BIOTECH PVT. LTD. An ISO 11134:2014, ISO 13485:2016, ISO 9001:2015 ← CIN NO.: U24210MH1995PTC095220S, S. No. 120/2, Laxmi Nagar, Umbarnala Road, Malkapur-443101, Dist.: Buldana (M.S.) India. Customer Care +91-8669083859 rdmsales@chaitanyagroupindia.com, mkt.cabt@chaitanyagroupindia.com, www.chaitanyagroupindia.com